**PUSH SPORTS ANKLE BRACE KICK**

The Push Sports Ankle Brace Kick stabilises the ankle joint and stops inversion (ankle sprain). With its innovative design and comfortable materials, the Push Sports Ankle Brace Kick offers the best possible protection against ankle sprain during sports activities.

**WHEN TO USE THE PUSH SPORTS ANKLE BRACE KICK?**

You have sprained your ankle before and want to avoid it happening again. In particular if you practise a sport with an increased risk of recurrent spraining, such as field or indoor sports. The Push Sports Ankle Brace Kick protects your ankle, significantly reducing the chances of new sprains.

**INDICATIONS**

- To prevent recurring inversion trauma
- In the event of ankle instability

**CONTRAINDICATIONS**

- There are no known disorders where the Push Sports Ankle Brace Kick should not be worn.

**CAUTIONS**

- Ensure that the Push Sports Ankle Brace Kick is not fitted too tightly. Use the brace only for the appropriate indications. Proper use of your Push Sports brace may require diagnosis by a doctor or therapist. If pain or complaints persist without a clear cause, it is recommended that you consult your GP, physiotherapist or sports physician.
- The Push Sports Ankle Brace Kick can be used in combination with an insole, provided that the space within the shoe allows such.
- For maximum comfort, it is recommended to wear the ankle brace with a sock.
- This product restricts inversion (inward tilting of the foot), but allows for largely unrestricted eversion (outward tilting of the foot).
- To make it easier to get into your shoe, it is recommended that you open up the laces of your shoe as far as possible.
- Push Sports braces were designed to be worn on an undamaged skin. Damaged skin must first be covered with suitable material.
- Optimal functioning can only be achieved by selecting the right size. Therefore, preferably fit the brace before purchasing, to check its size. In the case of an online purchase, always ensure that there is a good return service. Your reseller's professional staff will explain the best way to fit your Push Sports brace. We recommend that you keep this product information in a safe place, to be able to consult it for fitting instructions whenever needed.
- Check your Push Sports brace each time before use for signs of wear or ageing of components or seams. Optimal performance of your Push Sports brace can only be guaranteed if the product is fully intact.

- Within the framework of EU regulation 2017/745 on medical devices, the user must report any serious incident involving the product to the manufacturer and to the competent authority of the Member State in which the user is established.

**SPECIAL FEATURES**

The Push Sports Kick has been developed to restrict inversion (inward tilting of the foot) and allows for largely unrestricted eversion (outward tilting of the foot). In the event of inversion, the straps tighten and prevent the spraining movement. The brace's design allows for flexion and extension of the ankle joint, without unnecessary restriction during sports activities. There is minimal material on the instep of the foot, for unrestricted ball control. Its slender design makes the ankle brace eminently suitable for use inside sports footwear, including soccer cleats.

**TAKING MEASUREMENTS**

The Push Sports Ankle Brace Kick is available in three sizes, and in right and left versions. Select the right brace size on the basis of the circumference of your instep and your shoe size.

**FITTING**

- Open the three hook and loop fasteners and pull the brace over the foot. The heel section of the brace needs to fit snugly around the heel. Ensure the entire elastic base fits closely against the foot and is not twisted.
- Place the foot at a 90 degree angle and put the stiffener on the inside of the ankle in the correct position behind the malleolus (ankle bone) by tightening and fastening the small strap with the single blue mark. **Note:** The stiffener should not touch the malleolus.
- Then tighten and fasten the small strap with the two blue marks.
- Fasten the elastic strap to desired tightness around the ankle (Fig. 5).

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