




MD300CI216

FINGERTIP PULSE OXIMETER

ChoiceMMed Hospital Fingertip Pulse Oximeter Serial provides non-invasive, fast, accurate measurement for patients' oxygen saturation (SpO₂) and pulse rate (PR). Special product design makes it more durable and more suitable for hospital use. Real-time limit indication helps doctors and nurses catch the change of patients' vital signs on time.



FEATURES

- ♥ Colorful color OLED displays SpO₂, PR, Perfusion Index (PI)* and waveform.
- ♥ 2 display modes.
- ♥ 2pcs AAA-size alkaline batteries; battery-low indicator; real-time battery status indication.
- ♥ Automatically or manually power on, automatically power off.
- ♥ Weak or unstable signal prompt.
- ♥ Certification: CE, FDA



TECHNICAL SPECIFICATIONS

Application Range			
Adult, adolescent and child			
SpO ₂		Pulse Rate	
Measurement Range:	70%~100%	Measurement Range:	30bpm~250bpm
Resolution:	1%	Resolution:	1bpm
Measurement Accuracy:	70%~100% ±2%; 0%~69% no definition	Measurement Accuracy:	30bpm~99bpm ±2bpm; 100bpm~250bpm ±2%
PI*		Display	
Measurement Range:	0.1%~20%	Type:	OLED
Power Requirements		Dimensions and Weight	
Battery:	Two AAA-size 1.5V Batteries	Dimension:	61.5 (L) × 32 (W) × 31.5 (H) ±3mm
Power Consumption:	Less than 40mA	Weight:	60±5g (including batteries)
Environment Requirements			
Operation Temperature:	5°C~40°C	Operation Humidity:	15% ~ 93% no condensation
Storage Temperature:	-25°C~+70°C	Storage Humidity:	≤93%, no condensation

*PI: Perfusion Index, only for the devices sold in EU version.

Beijing Choice Electronic Tech Co., Ltd.

Room 4104, No. A12 Yuquan Road, Haidian District,
Beijing 100143, P.R.China
Website: www.choicemmed.com Fax: 86 10-88204632
Email: market@choicemmed.com Tel: 86 10-88204630/31

ChoiceMMed
Providing a Better Life

NOTE The specifications are subject to changes without prior notice.
Revised date: Jan 10, 2019