




MD300C55

FINGERTIP PULSE OXIMETER

ChoiceMMed Fingertip Pulse Oximeter for Children uses cartoon appearance design, reducing the fear of the children, and make them actively cooperate with medical staff checks, smaller size & less weight design will not cause too much pressure on kid's fingers, let Mom and Dad feel more ease.



FEATURES

- ♥ Cartoon appearance design
- ♥ Smaller size & less weight
- ♥ One button easy to operate
- ♥ Dual-color OLED screen displays SpO₂, PR, waveform, and pulse bar
- ♥ Adjustable level 1-10 brightness, 4-direction 6 display modes
- ♥ Unstable or weak signal prompt
- ♥ Real-time battery status indication
- ♥ 2pcs AAA-size batteries



TECHNICAL SPECIFICATIONS

Application Range			
adult, adolescent and child			
SpO ₂		Pulse Rate	
Measurement range:	70%~100%	Measurement range:	30bpm~250bpm
Resolution:	1%	Resolution:	1bpm
Measurement accuracy:	70%~100% ±2%; 0%~69% no definition	Measurement accuracy:	30bpm~99bpm ±2bpm; 100bpm~250bpm ±2%
Display		Dimensions and weight	
Type:	OLED	Dimension:	49mm (L) * 28mm (W) * 28mm (H)
Parameters:	SpO ₂ , PR, waveform, and pulse bar	Weight:	45g (including batteries)
Power Requirements			
Battery:	Two AAA-size 1.5V Batteries	Power consumption:	Less than 40mA
Environment Requirements			
Operation Temperature:	5°C~40°C	Operation Humidity:	15% ~ 93% no condensation
Storage Temperature:	-25°C~+70°C	Storage Humidity:	≤93%, no condensation
Accessories			
One silicone rubber boot	Two AAA batteries	One user manual	One lanyard One pouch

Beijing Choice Electronic Tech Co., Ltd.

Room 4104, No. A12 Yuquan Road, Haidian District,
Beijing 100143, P.R.China
Website: www.choicemmed.com Fax: 86 10-88204632
Email: market@choicemmed.com Tel: 86 10-88204630/31

ChoiceMMed
Providing a Better Life

NOTE The specifications are subject to changes without prior notice.
Revised Date: Jan 9, 2019